

MELTING POINT ATC PROBE MODEL No: ATC -11

Features:

- Class A simplex RTD sensor element.
- Three wire construction.



Specifications:

- a) Measurement Range : -50°C to 350 °C.
b) Element Type : RTD pt 100 Class A, (DIN 43760 Std.)
c) Accuracy : $\pm(0.15+0.002 * t)^\circ \text{C}$ Where 't' is the measured temperature.
d) Sensor Material : Glass
e) Sheath diameter : OD: 6mm
f) Cable Type : 7/34, AWG 3Core Teflon
g) Sensor Length : 170mm
h) Cable length : 1Mtr. Extn. Cable with male sockets.

Ordering Information:

S.No	Description	Code
1	ATC Probe (MP 96 Sensor)	ATC - 11

For Further Information Contact:

POLMON INSTRUMENTS PVT.LTD

(AN ISO 9001: 2000 Certified Company)

“POLMON HOUSE” Nizampet Road, Kukatpally, Hyderabad – 500 072

Ph No: 040-23053046, 23057308, 23055970, 23057642

Customer Care: 040-23892920 Fax: 040 – 23055425

e-mail: polmon@polmon.com URL: www.polmon.com