

CONDUCTIVITY CALIBRATOR

CAL-500



Conductivity Calibrator (Model : CAL-500) comes from POLMON ranges of Calibrators. The calibrator is designed to make calibrations convenient, quick and accurate.

Features

- Fixed precision resistors to eliminate drift
- 6 selectable Calibration Ranges
- Quick plug In for BNC & Screw terminal equipments
- ETDC certification (Optional)

Applications

- Conductivity instruments

Specifications

Technical		Physical	
Model No	CAL-500	Enclosure	ABS Plastic
Conductivity Ranges	100ms, 10ms, 1ms, 100µs, 10µs, 1µs	Dimensions	95X48X38mm
Resistivity Ranges	10Ω, 100Ω, 1K, 10K, 100K, 1M	Weight	0.135Kg
Accuracy	0.1% of the selected range	Environmental	
Selection	Manual	Operating Temperature	10 to 50 °C
Termination	BNC Clips	Humidity	10-95%RH, Non-Condensing

Ordering Information for CAL-500	
<i>Model No.</i>	<i>Range</i>
CAL500	01 – Standard (100 mS, 10 mS, 1mS, 100 µS, 10 µS, 1µS)